


Customer:  
 Customer Code:  
 Phone:

Contact:  
 Commission:

Lift chairs



| Pos. | Qty | Part Number | Designation                        |
|------|-----|-------------|------------------------------------|
| 1    | 1   | 69MECHASS1  | Frame                              |
| 2    | 1   | 68MRART2P1  | Articulated Frame (Left and right) |
| 3    | 1   | 68MRRELEV1  | Lift out System                    |
| 4    | 1   | 68MROTRL1   | Motor                              |
| 5    | 1   | 68MRTELRL1  | Remote control                     |
| 6    | 1   | 68MRTRALR1  | Transformer                        |

**Lift out Chair**



Lift out Chair

**Edition  
 06-02**